

## **NARKOMED AV2+ Ventilator**

### **TECHNICAL SERVICE BULLETIN # 171**

**Product:** Narkomed Anesthesia Systems with AV2+ Ventilators  
(Narkomed 2C, Narkomed 4, Narkomed GS, Narkomed Mobile,  
Narkomed Military)

**Status:** Class 3- Information

**Date:** July 23, 2003

**Re:** New Norgren Repair Kit

**Reference Doc:** ECN 03-0065

**Reason/Rationale:** To correct the occurrence of intermittent internal leakage across the regulator valve seat. This is observed as a continuous increase of pressure, as indicated on the Inspiratory Flow gauge, during expiration.

**Solution/Procedure:** Install new Norgren Repair Kit with stronger spring.

**Devices affected:** All AV2+ ventilators using Norgren Regulators (DMI Part number 4114252)

**Timing/When:** Upon next service related to Inspiratory Flow Regulator symptom described above under "Reason/Rationale".

**Additional Info:** This part **replaces** Norgren Repair kit part number 4115151.  
Part number 4115151 will no longer be available.

Narkomed MRI, MRI-2 Anesthesia Systems with AV2+ ventilators are excluded from using this repair part due to no serviceability below the top level assembly (to prevent damage to shielding).

**Cost:** Repair kit can be obtained at the then current list price.

**Ordering Info/Parts:** New Part number: **4118241- Kit- Repair- Norgren Reg**  
This can be order thru your Customer Care Area Representative.

**Distribution:** DrägerService, All Authorized Service Organizations

If you have any questions, please contact Technical Support by phone at 1-800-543-5047 or by fax at 1-215-723-5935.